CATTANI

air technology

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STERILIZING FILTER  
FILTRE BACTÉRIOLOGIQUE  
STERILISATIONSFILTER  
FILTRO ESTERILIZADOR

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


## **STERILIZING FILTER TREATMENT FOR COMPRESSED AIR**

The "SA" Balston filtering cartridge (B), with borosilicate glass microfiber filtering media, delivers a 99.9999+% efficiency on 0.01  $\mu$  m particles.

### **FUNCTIONING**

Aspirated air goes through the first paper filter (A) located on the head; after cooling, the air goes through the second filter (R) located at the base of the drying column; after the drying process through silica gel, there is the third filter in sintered bronze (H). At this stage the air enters the Balston cartridge (S) and moves on to the tank.

Balston filters can be sterilized in autoclave at 135  $^{\circ}$ C (max. 20 cycles); when used with clean air life of the filters is about 1 year, unless operating conditions require more frequent substitutions.

### **REPLACING STERILIZING FILTERS**

Unplug the unit, unscrew the two screws (C) on the column cover; use disposable gloves and place in autoclave immediately to avoid contamination.

### **REPLACING PAPER FILTERS WITH BALSTON FILTERS**

Unplug the unit, unscrew screws (C), remove OR (D), replace paper filter "F" with BALSTON filter.

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CATTANI S.p.A.